

How to install/uninstall/run PowerDicom


Installation

Note: For the installation of PowerDicom, administrative rights are necessary. Please contact your system administrator for assistance.

Run PowerDicom.msi




Click 'Next'




Select 'I accept', click 'Next'


How to Install/Uninstall/Run PowerDicom



Click 'Next'



PowerDicom will be installed.



Click 'Finish'

PowerDicom Installation part 1 is finished now.
You must register PowerDicom in order to use it. Proceed with chapter 'Registration' for the Basic, Professional, or Expert Edition or chapter 'Product Activation' for the Personal Edition respectively.

Registration

Please read from the following table which section applies to you:

License Type\ Edition	Personal	Basic	Professional	Expert
Single User	3. Product Activation	1. Registration with License File		
Concurrent	n.a.		2. Registration with License Server	

If you are using the **USB version** please make sure that the PowerDicom USB drive is mounted properly to your computer and proceed with chapter "Running PowerDicom".


If you are using a **Trial or Demo edition** of PowerDicom none of the above is applicable for you. Please proceed with chapter "Running PowerDicom".

1. Registration with License File

This section only applies to single user licenses of the Basic, Professional, and Expert Editions.

The Registration of PowerDicom will start automatically when running PowerDicom unregisteredly. For registration, there are no administrative rights necessary.

Registration can also be invoked from *Start -> program files -> PowerDicom -> Register now*




Order a PowerDicom License:

If you don't have a license, click 'Order License'.

If you have a valid license, please press 'Proceed' and proceed with chapter 'Finish Installation'

How to Install/Uninstall/Run PowerDicom

License file request dialog:




The image shows the 'PowerDicom Registrator' window. It has a title bar with the text 'PowerDicom Registrator'. Inside, there is a logo on the left and contact information on the right: Internet links to <http://www.MHGSoft.de/> and <http://www.dicom-solutions.com/>, and an email address [mailto: mhgs@MHGSoft.de](mailto:mhgs@MHGSoft.de). Below this is a section titled 'Your Registration Data:' containing several text input fields: 'Please type in your full name (first name, last name)', 'Please type in your eMail address', 'Please type in your company and department', 'Please type in your location/city', 'User name (account name)', and 'Computer name (host name)'. There are also radio buttons for 'Access Mode' (Read / Write (modifying files) and Read-only (not modifying files)), a checkbox for 'Single-Modality', a dropdown menu, and checkboxes for 'Expert Edition' and 'Demo'. At the bottom, there is a 'Field for Annotations:' with a large text area. A link 'Display privacy statement' is on the bottom left, and buttons 'Go', 'Help', 'About', and 'Exit' are on the bottom right.

Fill out the registration form. The eMail address field is for the very eMail address on which we will send you the license file. So please check its spelling thoroughly.

If you want to license for a single modality only, check the respective checkbox and choose a modality from the list. If the desired modality is not available, you may enter an expressive free text instead.


If you got a process number, please enter it in the 'Field for Annotations'.

Click on 'Go'



Click 'Prepare email'. Your personal registration email will be prepared.

How to Install/Uninstall/Run PowerDicom



Send this email to mhgs@MHGSoft.de to obtain a license.

You will receive your license file via email.

Save this license file to hard disk. You will need it for completing the installation process.


Please proceed with chapter 'Finish Installation'

Finish Installation

This section applies only to the single user licenses of the Basic, Professional, and Expert Edition.


You have a valid license file.

Run PowerDicom, e.g. *Start -> Program Files -> PowerDicom -> PowerDicom*



Click 'Proceed' to finish installation.

Next you must navigate to your license file.



Click 'Open'


Installation is finished now. PowerDicom is ready for use.

2. Registration with License Server

This section only applies to concurrent licenses of the Professional and Expert Editions.

Concurrent licenses require that a Licentaury License Server is available on your network. Please contact your network administrator for details.

The Registration of PowerDicom will start automatically when running PowerDicom for the first time. For registration, there are no administrative rights necessary.



Enter the name or the IP address of the computer, which hosts the Licentaury License Server, into the appropriate edit box. If you don't know the host name or IP address, please ask your administrator.

After clicking on OK PowerDicom will start and complete the registration process.


3. Product Activation

This section applies only to the **Personal Edition** of PowerDicom.

Product Activation will start automatically when running PowerDicom unregisteredly. For product activation, no administrative rights are required.

Select the License Data you received via eMail and copy them to the Clipboard (Ctrl+C).

Then paste (Ctrl+V) them to the Product Activation form and click on OK.



PowerDicom Personal Edition is now activated and ready for use.

Running PowerDicom

There are many ways to start PowerDicom:


- Start->Programs->PowerDicom->PowerDicom or
- Doubleclick PowerDicom Shortcut on your desktop or
- Doubleclick an Dicom File with extension *.ima, *.img, *.dcm or
- Select a Dicomfile, press right mousebutton and choose 'Send To' PowerDicom or
- Commandline: PowerDicom.exe *filename1* [*filename2 ...*]

Uninstall PowerDicom


Removing user specific files and settings.

Each user on a computer should remove changes in his profile prior to uninstalling PowerDicom from the computer:

Start RemovePD.exe, that is located in the PowerDicom directory:




Click on OK.



How to Install/Uninstall/Run PowerDicom


Removing PowerDicom from the computer

Run Uninstaller: Start -> Settings -> Control Panel -> Add or Remove programs



Select PowerDicom

Click 'Remove'



Click 'Yes'.

PowerDicom will be removed, all files, links, settings.